Hello.

QUICK START GUIDE

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CONTENTS

WELCOME	2
OVERVIEW	2
In The Box	2
Product Introduction	2
CAUTION	3
USING DR'S SECRET LIP	4
Application Modes	5
Notes	8
CLEANING	9
CHARGING	9
WARRANTY	11
SPECIFICATIONS	12

WELCOME

Congratulations on taking the first step towards a skin-renewing experience. Before you embark on your journey, please take a moment to read through the instructions in this User Guide.

OVERVIEW


DR's Secret LIF is a revolutionary device that when used as part of your daily beauty regime, helps to cleanse your skin and prep it for your skin care ritual.

IN THE BOX

This box contains:

1. DR's Secret LIF (Main Unit)
2. Power Adaptor
3. USB Cable
4. User Guide

PRODUCT INTRODUCTION

CAUTION

- The device is only for its intended use as described in this User Guide.
- Use only with dry hands. DO NOT use while bathing or in the shower. DO NOT rinse device with water or submerge in water.
- Remove all jewelry such as rings, necklaces and bracelets before use.
- DO NOT use if you have skin infection or malignancy or sunburn.
- DO NOT use if you are on steroid treatment or feeling feverish.
- If you are pregnant or suspect that you are pregnant, please consult your physician before use.
- Please note that LIF emits low level electric currents during use. DO NOT use LIF if you are fitted with a pacemaker or similar device.


If you have cardiac or heart conditions, or have a metal implant, metallic dental fixture or filling in your body, please consult your doctor before use.

- Keep out of reach of children. DO NOT use on children below 12 years of age.
- Avoid leaving your device in direct sunlight and DO NOT expose to extreme heat or in humid areas.
- Keep the device away from magnetic objects such as hard disks, video tapes, computer monitors and televisions.
- Never use the Power Adaptor or USB Cable if damaged. Make sure the charging accessories and device is dry when charging. DO NOT charge device near water.
- Handle the product with care. DO NOT drop the device.

USING DR'S SECRET LIF

The 6 self-timed modes:

- C** **Cleansing Mode** ~ ~ +/—
- V** **Vitamin C Mode** ~ ~ +/—
- B** **Brightening Mode** ~ ~ +/— ~
- L** **Lifting Mode** ~ ~ +/— ~
- ~ **Heat Mode** ~ ~ ~
- ~ **Light Mode** ~ ~ ~

It is recommended to use **B** or **L Modes** daily. Use the full regime **C-V-B-L** 2 to 3 times weekly or whenever you use Vitamin C serum. Use your Vitamin C serum with **V Mode**, apply your other skin care products and end with **B** or **L Modes**.

APPLICATION MODES

- During use, please always make sure there is sufficient skin care product applied evenly on the skin to allow smooth gliding of DR's Secret LIF over the face.
- For modes with ~ **Heat**, DO NOT pause and hold at one spot for too long. Pause and hold for only 1-2 seconds. You may also use gliding motion.
- It is normal for skin to feel warm or look slightly reddened after using modes with ~ **Heat**. You can apply ~ **Blue Light** or any soothing skin care product of your choice if needed.

C Cleansing Mode

3 minutes, composed of:

- 1.5 minutes : ~ ~ +/—
- 1.5 minutes : ~ ~ +/—

Cleansing Mode helps to pull out impurities, dirt and remove residual makeup left on skin.

Complementary products: Aestier Cleansing Cream A1.

V Vitamin C Mode

3 minutes: ~ ~ +/—

Vitamin C Mode helps your skin to absorb the antioxidant more efficiently. Vitamin C can also promote skin renewal.

Complementary products: DR's Secret Toner T2 followed by DR's Secret C15 Essence 7.

B Brightening Mode

7 minutes, composed of:

- 3 minutes : ~ ~ ~
- 2 minutes : ~ ~ +/—
- 2 minutes : ~ ~ +/—

Brightening Mode helps to promote skin radiance and improve skin texture. It also prepares your skin for absorption of your skin care product. You may also use Brightening Mode for eye treatment.

Complementary products: DR's Secret Skinlight T3, DR's Secret Skinrecon T4, Aestier Emulsion A6/DR's Secret Aqua Boost Serum 10 or DR's Secret Toner T2 followed by DR's Secret C15 Essence 7. For eye treatment, use Aestier Eye Cream.

L Lifting Mode

5 minutes, composed of:

- 3 minutes : ~ ~ ~
- 2 minutes : ~ ~ +/—

Lifting mode helps to sculpt your facial contours for a defined appearance. You may also use Lifting Mode for eye treatment.

Complementary products: Aestier Emulsion A6/DR's Secret Aqua Boost Serum 10 or DR's Secret Toner T2 followed by DR's Secret C15 Essence 7. For eye treatment, use Aestier Eye Cream.

🌊 Heat Mode

15 minutes: 🌊 🌊

Heat Mode helps to target sagging lines and fine lines. It can be used to target specific skin problems.

Complementary products: Aestier Emulsion A6/DR's Secret Aqua Boost Serum 10 or DR's Secret Toner T2 followed by DR's Secret C15 Essence 7.

💡 Light Mode


15 minutes: 💡 💡

Light Mode allows users to enjoy convenient light therapy. It can be used to target specific skin problems.


Complementary products: DR's Secret Skinlight T3, DR's Secret Skinrecon T4, Aestier Emulsion A6/ DR's Secret Aqua Boost Serum 10, DR's Secret Spot Serum 8 or DR's Secret Refining Serum 9.

*Our products may be updated from time to time. Product information shown here may vary from the product packaging depending on time and region of purchase.

1. For **C Mode**, moisten face and apply cleansing cream to skin.
For **V Mode**, apply Vitamin C serum evenly on moist face.
For **B, L, 🌊 Heat** or **💡 Light Modes**, apply facial product evenly over face.
For eye treatment, apply eye cream evenly over eye area.
2. With dry hands, switch on DR's Secret LIF and choose the appropriate mode.
For **💡 Light Mode**, press once for **Red Light**, twice for **Yellow Light** and thrice for **Blue Light**.
3. Choose your intensity level. For eye treatment, please use only level 1 or 2.
4. **Glide** DR's Secret LIF over your face as shown.


For eye treatment, you may **glide** DR's Secret LIF over eye area or use a **Pause-and-Hold** method as shown. Please note to only pause and hold DR's Secret LIF on each spot for **only 1-2 seconds**.

For **💡 Light Mode**, you may **glide** DR's Secret LIF over your problem areas or use a **Pause-and-Hold** method as shown.

* Use only **Pause-and-Hold** method for DR's Secret Spot Serum and Refining Serum. *

5. Press the appropriate mode again to pause the program. Press again to resume.
For **Light Mode**, there is no pause function.
6. Apply more facial product or cleansing cream if needed.
7. DR's Secret LIF will beep to indicate completion of mode. Light on device head will be switched off.
8. Tissue off gently or rinse with water (for **C Mode** only).
9. Wipe the device head with a wet towel or alcohol-free wet tissue. Dry with cotton or tissue. Pay extra attention to wipe the grooves on the device head.

* For cystic acne use **Blue Light** until acne clears up. *

NOTES

- For first-time users, please start from level 1 and increase slowly. It is normal to experience mild warming sensations during **B, L** and **Heat Modes**.
- If you feel any discomfort, please lower the intensity level, or apply more serum on skin. If discomfort persists, please consult your dermatologist.


CLEANING

The life span of DR's Secret LIF will be greatly increased with proper care and maintenance.


- Ensure that the device is detached from all electrical outlets before cleaning the device head.
- Wipe the device head with a wet towel or alcohol-free wet tissue after each use. Dry with cotton or tissue. Pay extra attention to wipe the grooves on the device head.
- DO NOT rinse your device under water or submerge into water.

CHARGING

- DR's Secret LIF takes 1.5 hours to be fully charged. Please charge the device before using.
- There are 2 methods to charge DR's Secret LIF.
 - Connect USB Cable to device


- Option 1: Connect to Power Adaptor and charge

- Option 2: Connect to USB port (e.g., to computer)

* Please note that charging time will take up to 4.5 hours when charging via computer. *

- The battery indicator light shows the battery status.
 - Blinking red light: Battery is close to being fully discharged. Device may automatically shut down during use.
 - Red light: Battery level is low. Please charge device when possible.
 - Yellow light: Battery level is moderate.
 - Green light: Battery level is full.
- When the battery is fully discharged, the device will beep continuously for 2 seconds. You will not be able to switch on the device.
- During charging, the battery indicator light will automatically light up to show the battery status.
 - Blinking red light: Battery is fully discharged. Device will not be able to switch on.
 - Red light: Battery level is low.
 - Yellow light: Battery level is moderate.
 - Green light: Battery level is full.
- You may also continue to use DR's Secret LIF while it is charging. However, please remember not to hold the device with wet hands while it is charging.
- You do not need to charge DR's Secret LIF continuously.

WARRANTY

- Your DR's Secret LIF is fully guaranteed against manufacturing defects for a period of 12 months from the date of purchase unless stated otherwise.
- Local warranty applies to the country of the original purchase only.
- Original tax invoice must be presented as proof of purchase to our service staff when DR's Secret LIF is sent for servicing/repair.
- This warranty is not transferable.
- This warranty can also be activated online by following the steps below:
 1. Go to 'www.bwlgroup.com'
 2. Choose your location
 3. Scroll down and click on 'Device Product Warranty' under 'Quick Links' in the website footer
 4. Enter your details and product individual serial number
 5. Click 'Submit'
- * Best World Lifestyle recognises online warranty activation date as date of purchase, up to a maximum difference of 6 months from tax invoice date. When online warranty activation date exceeds 6 months from tax invoice date, a maximum of 6 months is added to tax invoice date. That date will be taken as date of purchase. *
- This warranty shall immediately cease and made void if the product is found to have been modified, tampered with or repaired by an unauthorised person.
- This warranty shall be subjected to the following limitations and exclusions:
 - Malfunctions or damages resulting from natural disasters, fire, civil unrest and/or accidents
 - Defects, misuse or damage by negligence and abuse
 - Wear and tear of accessories

- Damage during transportation to/from BWL retail centre
- Product has not been maintained or operated in accordance with the instructions stated inside User Guide
- Best World Lifestyle reserves the right to amend the above terms and conditions without prior notice.

SPECIFICATIONS

Power source	100~240VAC / 50~60Hz
Power consumption (main unit)	7.2W MAX
Power output (adaptor)	DC5V, 2.0A
Supply currency	Modified Square / Sine Wave
Battery	Li-ion 3.7V, 1330mA/h
Dimension and weight (main unit)	Dimensions: 55mm*165mm*43mm Weight: 149g
Electrical intensity at lowest level	43V rms
Electrical intensity at highest level	61V rms

Disclaimer

While every effort has been made to ensure that information contained in this User Guide is accurate, BWL shall not be liable for any errors therein.

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DR's Secret and BWL are registered trademarks of Best World.
© 2022 Best World. All rights reserved.